


# TEG Manager® Analysis


Patient ID:	Test Start Time: 3/10/2022 1:59 PM Duration: 52 min
Segment:	Process:
Test Information:	Test Status: Test Completed Operator ID: SAMIRRUDDING
	Cartridge: CM Citrated K,KH,RT,FF Printed Cartridge Lot: 072020-1
	Device Name: TEG 6s Device Location: <span style="background-color: black; color: black;">XXXXXXXXXX</span> Device SN: T1-20G-102075



	TEG-ACT (sec)	R (min)	K (min)	ANGLE (deg)	A10 (mm)	MA (mm)	LY30 (%)
<b>CK</b>		6.0 4.6 - 9.1	1.0 0.8 - 2.1	75.9 63 - 78		68.2 52 - 69	0.0 0.0 - 2.6
<b>CRT</b>	69.2 82 - 152	0.2 0.3 - 1.1	0.8 0.8 - 2.7	80.2 60 - 78	66.4 44 - 67	68.3 52 - 70	0.0 0.0 - 2.2
<b>CKH</b>		6.1 4.3 - 8.3	1.0 0.8 - 1.9	76.6 64 - 77		68.2 52 - 69	
<b>CFF</b>					34.6 15 - 30	36.9 15 - 32	

# TEG Manager® Analysis

Patient ID:	Test Start Time: 3/14/2022 1:43 PM Duration: 54 min
Segment:	Process:
Test Information:	Test Status: Test Completed Operator ID: JAROMIN
	Cartridge: CM Citrated K,KH,RT,FF Printed Cartridge Lot: 111020-1
	Device Name: TEG 6s Device Location: <span style="background-color: black; color: black;">XXXXXXXXXX</span> Device SN: T1-20G-102075



	TEG-ACT (sec)	R (min)	K (min)	ANGLE (deg)	A10 (mm)	MA (mm)	LY30 (%)
<b>CK</b>		5.2 4.6 - 9.1	1.1 0.8 - 2.1	74.9 63 - 78		64.0 52 - 69	0.0 0.0 - 2.6
<b>CRT</b>	78.5 82 - 152	0.3 0.3 - 1.1	0.8 0.8 - 2.7	79.2 60 - 78	63.4 44 - 67	65.8 52 - 70	0.0 0.0 - 2.2
<b>CKH</b>		5.2 4.3 - 8.3	1.0 0.8 - 1.9	76.0 64 - 77		65.0 52 - 69	
<b>CFF</b>					27.5 15 - 30	28.9 15 - 32	